PROUD SPONSOR

WE NEED THIS TO BE

## Easy to use

...BUT WE DON'T WANT TO GIVE UP ANY

## Performance

NATHAN YATES UNIVERSITY OF PITTSBURGH

- Plug & Spray operation
- Single connection
- Temperature-controlled
- Factory-verified voltage connections
- Chip-to-chip consistency
- Inexpensive to operate



PicoChip Stage
Two PicoChip Columns

\$99900

You asked for it. We listened.




NANOSPRAY COLUMNS

The market-proven high-performance PicoFrit® column, with a factory-configured high-voltage liquid-junction connection integrated into an easy to use, pre-optimized format. The PicoChip is not only easy to use, but easy on the budget and it's all backed by New Objective's commitment to service and quality.

Model PCH-550 For Thermo LTQ




Place PicoChip on source




Tip is protected inside chip



Insert chip key to expose tip



Tip is exposed to spray



Collect data



New Objective, Inc.
Two Constitution Way
Woburn, MA 01801-1023

888 220 2998 Orders 781 933 9560 Phone 781 933 9564 Fax www.newobjective.com sales@newobjective.com

Innovation in High-Sensitivity LC-MS Research

